

# The Twinkle Factory

## match<sub>550</sub>

a.k.a. HBR-2,5DM

Reference match550-50X

Quantity 66 µg

Store at 2-8 °C


### Spectral properties of match<sub>550</sub> when FIRE<sub>tag</sub> & FIRE<sub>mate</sub> interact

Excitation wavelength 500 nm

Emission wavelength 550 nm

**match<sub>550</sub>** is a fluorogenic molecular glue that can be used to selectively induce the dimerization of proteins fused to cognate FIRE<sub>mate</sub> and FIRE<sub>tag</sub>. One vial includes 250 nmol of **match<sub>550</sub>**, enabling to prepare 50 mL of a 5 µM dimerizing solution.


The Twinkle Factory technology CATCHFIRE (Chemically Assisted Tethering of Chimera by Fluorogenic Induced Recognition) is a fluorogenic chemically induced dimerization technology that enables one not only to artificially control the proximity of two proteins of interest in cells, but also to see their interactions. The two proteins are genetically fused to two small protein domains (FIRE<sub>mate</sub> et FIRE<sub>tag</sub>) which are capable of interacting together in a reversible fashion in the presence of small fluorogenic molecular glues such as **match<sub>550</sub>**. When the two domains interact, the fluorescence of **match<sub>550</sub>** increases 100 fold, enabling to see the newly induced interaction by fluorescence microscopy.



The use of CATCHFIRE implies cloning and expressing proteins fused to FIRE<sub>mate</sub> and FIRE<sub>tag</sub>, and treating cells with **match<sub>550</sub>**. The protocol is described below. Note that proteins of interest can be expressed with FIRE<sub>mate</sub> and FIRE<sub>tag</sub> as either N- or C-terminal fusions.

The Twinkle Factory provides a range of fluorogenic dimerizers of various emission wavelengths, **match<sub>540</sub>**, **match<sub>550</sub>**, & **match<sub>600</sub>**, and a non-fluorogenic dimerizer, **match<sub>dark</sub>**. All of them are compatible with FIRE<sub>mate</sub> and FIRE<sub>tag</sub>.

Cells expressing proteins fused to FIRE<sub>mate</sub> and FIRE<sub>tag</sub> are not supplied by The Twinkle Factory. Plasmids containing FIRE<sub>mate</sub> and FIRE<sub>tag</sub> genes are available at Addgene [www.addgene.org/Arnaud\\_Gautier/](http://www.addgene.org/Arnaud_Gautier/).



# The Twinkle Factory

## Protocol of labeling in living cells

Dissolve one vial of **match<sub>550</sub>** in 50  $\mu\text{L}$  of DMSO to yield a 5 mM stock solution. Mix by vortexing for few seconds until all the **match<sub>550</sub>** is dissolved. Note that different stock concentrations can be made depending on your requirements. **match<sub>550</sub>** is soluble in DMSO up to at least 50 mM.

Dilute the stock solution 1:1000 in medium or buffer to yield a 5  $\mu\text{M}$  dimerizing solution. Mix thoroughly by vortexing. For best performance, add **match<sub>550</sub>** to serum-free medium or buffer, and do not keep/store the dimerizing solution. Note that different concentrations can be made depending on your requirements. Optimal concentrations range from 1 to 10  $\mu\text{M}$ .

Remove the cell culture medium, wash with D-PBS, and replace the buffer with the dimerizing solution to induce protein dimerization. Image with appropriate settings.

To reverse the dimerization, remove the dimerizing solution, wash with D-PBS, and replace with culture medium.

## Storage

Dry **match<sub>550</sub>** should be stored at 2-8  $^{\circ}\text{C}$  in the dark. Once dissolved in DMSO, the solution should be aliquoted to avoid repeated freeze/thaw cycles and stored at  $-20^{\circ}\text{C}$  in the dark. With proper storage, **match<sub>550</sub>** should be stable at least three years dry or 6 months dissolved in DMSO.

## Purity and Characterization

Purity of **match<sub>550</sub>** was determined to be  $> 99\%$  by nuclear magnetic resonance (NMR) and elementary analysis.

## References

A fluorogenic chemically induced dimerization technology for controlling, imaging and sensing protein proximity. *Nature Methods* **20**, 1553–1562 (2023). [doi.org/10.1038/s41592-023-01988-8](https://doi.org/10.1038/s41592-023-01988-8)

## Notice to Buyer/User

The Buyer/User has a non-exclusive license to use this system or any component thereof for research use only. The products and/or their use may be covered by one or more of the following patents and patent applications:

- EP 3,164,411; JP 2017-527,261; US 10,138,278 (Fluorogen activating and shifting tag (FAST))
- EP 3,719,007; US 2022-0,169,682 (Split photoactive yellow protein complementation system and uses thereof)
- EP Appl. 22,306,308.2 (Proximity inducing system and uses thereof)

These patents and patent applications are owned by CNRS (France), PSL Research University (France), Sorbonne Université (France), and/or Institut Curie (France), and exclusively licensed to Twinkle Bioscience S.A.S., 45 rue d'Ulm, 75005 Paris, France. Commercial use of this system or any components thereof requires a license. For more information [contact@the-twinkle-factory.com](mailto:contact@the-twinkle-factory.com).