



The Twinkle Factory

Stain different, tag FAST.

^{TF}Carmine

a.k.a HPAR-3,5DOM

Reference 636715-250

Quantity 250 nmol

Store at 2-8 °C


Properties of ^{TF}Carmine when bound to nirFAST

Excitation wavelength	636 nm
Emission wavelength	715 nm
Molar absorption coefficient	41 mM ⁻¹ cm ⁻¹
Fluorescence quantum yield	7 %
Affinity constant at 25° C	0.036 μM

^{TF}Carmine is a membrane-permeant fluorogenic ligand that can be used to selectively label nirFAST-tagged proteins in solution, in living cells and in fixed cells. ^{TF}Carmine is almost non-fluorescent when free in solution, but strongly fluoresces when bound to nirFAST. ^{TF}Carmine was designed together with a variant of FAST, nirFAST, specifically for near-infrared labeling. It exclusively works with this variant and should not be used with FAST1 nor FAST2, pFAST, frFAST. Also, it is not recommended to split nirFAST for protein-protein interaction reporting while the affinity of nirFAST for its cognate ^{TF}Carmine is very high. This package includes 250 nmol of ^{TF}Carmine, enabling to prepare 50 mL of a 5 μM labeling solution.

The Twinkle Factory labeling technology enables the specific fluorescent labeling of any protein of interest. It is based on the instantaneous formation of a fluorescent molecular assembly between the small (14 kDa) protein tag FAST and fluorogenic ligands, ^{TF}Fluorogens. ^{TF}Fluorogens strongly fluoresce only when bound to FAST, enabling to detect and image FAST-tagged proteins with high contrast without the need of washing the excess of fluorogenic ligands. The labeling of FAST-tagged proteins with a ^{TF}Fluorogen is non-covalent and can be reversed if necessary by washing.

The Twinkle Factory labeling technology implies cloning and expressing of the FAST-tagged protein, and labeling the resulting fusion with the ^{TF}Fluorogen of choice. The labeling of FAST-tagged proteins is described below. Cells expressing FAST-tagged proteins are not supplied by The Twinkle Factory. Note that proteins of interest can be expressed with FAST as either an N- or a C-terminal fusion.



Absorbance (dotted line) and emission (solid line) spectra of ^{TF}Carmine bound to nirFAST

Protocol of labeling in living cells

Dissolve one vial of ^{TF}Carmine in 50 μL of DMSO to yield a 5 mM stock solution. Mix by vortexing for few seconds until



all the ^{TF}Carmine is dissolved. Note that different stock concentrations can be made depending on your requirements. ^{TF}Carmine is soluble in DMSO up to at least 50 mM.

Dilute the stock solution 1:500 in medium or buffer to yield a 10 μ M labeling solution. Mix thoroughly by vortexing. For best performance, add ^{TF}Carmine to serum-free medium or buffer, and do not keep/store the labeling solution. Note that different concentrations can be made depending on your requirements. Optimal concentrations range from 1 to 10 μ M.

Remove the cell culture medium, wash with D-PBS, and replace the buffer with the labeling solution. Incubate for 15-30 seconds and image the cells directly.

Image the cells using appropriate settings. nirFAST-tagged proteins labeled with ^{TF}Carmine have an excitation maximum at 636 nm and an emission maximum at 715 nm.

To reverse the labeling, remove the labeling solution, wash with D-PBS, and replace with culture medium. Or you can add ^{TF}Darth to the culture medium.

Protocol for labeling in fixed cells

Cells expressing nirFAST-tagged proteins can be fixed before labeling with standard fixation methods such as paraformaldehyde, ethanol, methanol. Once the fixation is performed, wash cells with D-PBS, and replace the buffer with a labeling solution (prepared in D-PBS). Incubate for 15-30 seconds and image the cells directly as above. To reverse the labeling, remove the labeling solution and wash with D-PBS.

Storage

Dry ^{TF}Carmine should be stored at 2-8 °C in the dark. Once dissolved in DMSO, the solution should be aliquoted to avoid repeated freeze/thaw cycles and stored at – 20 °C in the dark. With proper storage, ^{TF}Carmine should be stable at least three years dry or 6 months dissolved in DMSO.

Purity and Characterization

Purity of ^{TF}Carmine was determined to be >99% by nuclear magnetic resonance (NMR) and elementary analysis.

References

bioRxiv 2024-04 (doi.org/10.1101/2024.04.05.588310)

Notice to Buyer/User

The Buyer/User has a non-exclusive license to use this system or any component thereof for research use only. Commercial use of this system or any components thereof requires a license from Twinkle Bioscience S.A.S, 45 rue d'Ulm, 75005 Paris, France. For detailed information, e-mail contact@the-twinkle-factory.com. The products and/or their use may be covered by one or more of the following patents and patent applications:

- EP 3 164 411; JP 2017-527,261; US 10,138,278 (Fluorogen activating and shifting tag (FAST))

These patents and patent applications are owned by CNRS (France), PSL Research University (France), Sorbonne Université (France), and/or Institut Curie (France), and exclusively licensed to Twinkle Bioscience.